



Easy Pour Breast Milk Storage Bags

Key Features & How to Use



Patented design

Wide Pour in (1) and out (2), spouts for an easy transfer

10x better protection of the breast milk*

*through better oxygen barrier vs Medela's former bag

Flat freezing

Large capacity: 210 ml / 7 oz

Self-standing

Anti-spill anti-mess



Easy to write on label

Inverted for a more comfortable writing experience

Sanitized

Sterilized with e-beam treatment for hygiene, safely guarded in the transport pouch and ready to use

Recyclable product

Without rinsing, for extra peace of mind

Recyclable retail packaging

Made of cardboard

Recyclable transport pouch



No leaks

Patented design : ensure no leaks, self standing by design to avoid any spilling



x10 better protection* of breast milk

*Through a better oxygen barrier vs former's Medela bag

Double layer / Ultra Thick Material for an optimum protection.



Flat-freezing

Space optimization in the freezer or in the fridge.





Bigger capacity for a flatter freezing : 210 ml vs 180 ml



Recyclable and sanitized

Ready to use product, sterilized with e-beam treatment for hygiene

Easy to recycle : no need to wash or rinse the product before throwing away.



Storage

Storage place	Room temperature 16 °C to 25 °C	Refrigerator 4 °C or colder	Freezer -18 °C or colder	Previously frozen breast milk thawed in the refrigerator
Safe storage time	Up to four hours is best	Up to three days is best	Up to six months is best	Up to two hours at room temperature